

Salivary Steroid ELISA Kits

FDA approved

Cortisol

- Sensitivity: 1.14 ng/ml
- Standard range: 2 - 80 ng/ml

FDA approved

17 α -OH Progesterone

- Sensitivity: 0.003 ng/ml
- Standard range: 0.005 - 1 ng/ml

Androstendione

- Sensitivity: 5 pg/ml
- Standard range: 20 - 1000 pg/ml

FDA approved

Progesterone

- Sensitivity: 0.0055 ng/ml
- Standard range: 0.01 - 5 ng/ml

FDA approved

Testosterone

- Sensitivity: 0.006 ng/ml
- Standard range: 0.01 - 5 ng/ml

Free Estriol

- Sensitivity: 2 pg/ml
- Standard range: 0.1 - 12 pg/ml

FDA approved

Estradiol

- Sensitivity: 0.05 pg/ml
- Standard range: 0.1 - 100 pg/ml

Estrone*

- High Sensitivity
- Short incubation time

Cortisol HS (high sensitive)

- Sensitivity: 0.015 ng/ml
- Standard range: 0.1-30 ng/ml

DHEA

- High Sensitivity: 2.186 pg/ml
- Standard range: 10 - 1440 pg/ml

* Coming soon


Cardiac Diseases

PAPP-A (ultrasensitive)*

Pregnancy-Associated Plasma Protein A, a new candidate marker for the early diagnosis of **Acute Coronary Syndroms (ACS)**.

- Incubation time: 2 h 45 min.
- Standard range: 0 - 450 ng/ml
- Sample: 100 μ l Serum

*) Not available in USA



Liver Fibrosis

CTGF*

(Connective Tissue Growth Factor)
Potential new marker for liver fibrosis in early stages (F1, F2)

- Specificity 78% (liver fibrosis F1)
- Sensitivity 83,7% (liver fibrosis F1)

*) coming soon